



# Trilogy® Acetabular System

Surgical Technique



Versatility in a proven design





**Trilogy Acetabular System  
Surgical Technique****Table of Contents**

|   |          |
|---|----------|
| <b>Acetabular Reaming</b>                     | <b>2</b> |
| <b>Component Sizing</b>                       | <b>2</b> |
| <b>Positioner Assembly</b>                    | <b>3</b> |
| <b>Positioning and Alignment</b>              | <b>3</b> |
| <b>Implant Orientation</b>                    | <b>4</b> |
| <b>Implant Insertion</b>                      | <b>4</b> |
| <b>Screw Insertion</b>                        | <b>5</b> |
| <b>Liner Positioning</b>                      | <b>5</b> |
| <b>Liner Insertion</b>                        | <b>6</b> |
| <b>Locking Ring Replacement (if required)</b> | <b>6</b> |
| <b>Final Reduction</b>                        | <b>7</b> |
| <b>Ordering Information</b>                   | <b>8</b> |

## Acetabular Reaming

Use reamer to prepare acetabular implant site. Hold reamer steady and apply pressure in the same direction that the prosthesis will be implanted.



## Component Sizing


Surgical judgment is used to assess proper fit for press-fitting to avoid fracturing the acetabulum. Assess bone quality and determine appropriate implant size and type. Place Provisional Shell into the prepared acetabulum using the positioner to check the sizing. Do not impact the provisional. The provisional shell should protrude above the periphery of the acetabulum. Another method is to measure the “gap” through the polar hole of the provisional at the dome. The Provisional Liner can be inserted to perform a trial reduction.

**Note:** Careful assessment of this step is necessary, especially when using the non-holed, uni-holed, or spiked cup.




**Warning:** The *Trilogy* Acetabular Cup is a full hemisphere, and adequate bone stock is necessary for a press-fit application.

**Note:** A cup sized 2mm over the reamed preparation (size of last reamer shell used) will provide a 2mm press-fit.



Moderate protrusion above the periphery



Provisional Liner and Shell

**Warning:** Do not impact onto the provisional liner as damage to the liner may occur.

Optional MIS Offset Reamer Handle  
00-7804-080-00




## Positioner Assembly

Thread the implant onto the positioner securely. Position alignment connector as shown. Alignment connector should be in line with locking ring access. Impact alignment connector to engage onto the taper. Tighten Locking Nut. Impaction will be applied primarily to the dome of the implant.


*Alternatively: The Rim Impactor can be utilized. Remove the cap from the Cup Positioner and replace with the Rim Impactor. The Rim Impactor size should correspond with the implant size selected. The implant should be threaded onto the positioner securely. Impaction will be applied primarily to the rim of the implant.*

**Warning:** Use caution to avoid impacting the anti-rotational tabs on the rim of the implant.



## Positioning and Alignment

Ensure the patient is in the correct position. With the final prosthesis connected to Acetabular Cup Impactor-Positioner, insert shell into prepared site. Keep slap hammer portion of positioner vertical, and "A-Frame" or "Gunsight" extension parallel with patient's body to achieve 45 degrees abduction and 20 degrees forward flexion.






## Implant Orientation


Neurologic and vascular injury can be minimized if the posterior quadrants are used for trans-acetabular screw placement.\* The cluster-holed shell (A) should be positioned to allow screw placement in the posterior superior and/or posterior inferior quadrants of the acetabulum. If the spiked shell (B) is selected, the general position of the locking ring will be 10 o'clock or 2 o'clock. To ensure avoidance of the anterior inferior acetabulum, place the anterior spike in the pubic area.

\* Wasielewski RC, Cooperstein LA, Kruger MP, Rubash HE. Acetabular anatomy and the transacetabular fixation of screws in total hip arthroplasty. *J Bone Joint Surg.* 1990;72-A,4,501-508.



## Implant Insertion

The “Gunsight” or “A-Frame” Alignment Guides can be used to aid in cup positioning. Use short strikes with a small mallet to impact cup. The impact required to safely seat the implant is dictated by the bone quality. Due to impaction, it may be necessary to periodically re-tighten the collet or locking nut. Remove “A-Frame” or “Gunsight” device during continued impaction. Use screw holes, or the positioner hole on the spiked shell, to determine if the shell is fully seated.






Optional Offset Shell Inserters  
00-7804-025-20

## Screw Insertion


If screw placement is desired, drill a pilot hole by placing the drill through the Drill Guide in the desired screw hole (A). Measure the hole's depth with the Depth Gauge (B). Select the appropriate length screw and insert it in the hole with the Screwdriver (C). Place additional screws as necessary.

**Warning:** Avoid screw placement through the shell into the anterior inferior quadrant of the acetabulum to prevent injury to the neuro-vascular structures.




## Liner Positioning


Insert provisional liner into prosthesis to determine appropriate liner elevation. The anti-rotational slots on the trial liner should line up with the two anti-rotational tabs and black etch marks located on the rim of the shell. **Clear all soft tissue** from around the perimeter of the shell and assess visualization. It may be necessary to adjust retractors to optimize easy insertion of the prosthetic liner. To achieve complete congruency, dry inner shell thoroughly. Insert the final poly liner with liner positioner, or by hand. Align with anti-rotational tabs. Ensure that both tabs are engaged with prosthesis liner slots.



Anti-rotational tab and slot misaligned



Correct vertical alignment




Shell and polar boss misaligned


## Liner Insertion

Push the liner into the shell. Impaction may be necessary to fully engage liner into shell. Watch for the locking tabs to open and then close together. There will be a 2mm to 3mm gap between the locking ring tabs. The prosthetic liner is fully seated when the locking ring tabs can be moved from side-to-side with a probe. The ring serves to lock the liner in place. It must float freely inside the shell to function properly.

A  
Before Insertion




B  
Insertion



C  
Locked





D  
Fully Seated (side-to-side check)



## Locking Ring Replacement (if required)

If upon inspection, it is determined that the locking ring is not functioning properly or has become damaged, it must be replaced. If seated, the liner is removed using the Disassembly Device to open the ring, while at the same time using the Liner Extractor to clamp and lift the liner out of the shell.

To remove the old ring, push the ring to the side of the shell and grab the exposed portion on the opposite side of the shell. A firm and steady grip is required. To insert a new ring, first insert the tabs of the ring into the locking ring window. Then, press the ring into place. When properly seated, the ring will first pop into place and then float freely in the shell. Liner assembly can now take place.





## Final Reduction

Perform final reduction to assess range of motion, hip stability, and limb length. If “fine tuning” or a different liner style is needed, use the Disassembly Device to remove the liner.



## Ordering Information

### Trilogy Acetabular Shells

SHELLS are available in 2mm OD increments  
6200-SHELL OD-SHELL TYPE

| Prod. No.            | Description                               |
|----------------------|---|
| <b>Spiked</b>        |   |
| 00-6200-040-23       | F/M Acetabular Shell, 40mm OD, Spiked     |
| Through ↓            | Through ↓                                 |
| 00-6200-070-23       | F/M Acetabular Shell, 70mm OD, Spiked     |
| 65-6200-040-23       | HA/TCP Acetabular Shell, 40mm OD, Spiked  |
| Through ↓            | Through ↓                                 |
| 65-6200-070-23       | HA/TCP Acetabular Shell, 70mm OD, Spiked  |
| <b>Non-Holed</b>     |   |
| 00-6200-040-21       | F/M Acetabular Shell, 40mm OD, Solid      |
| Through ↓            | Through ↓                                 |
| 00-6200-070-21       | F/M Acetabular Shell, 70mm OD, Solid      |
| 65-6200-038-21       | HA/TCP Acetabular Shell, 38mm OD, Solid   |
| Through ↓            | Through ↓                                 |
| 65-6200-070-21       | HA/TCP Acetabular Shell, 70mm OD, Solid   |
| <b>Cluster-Holed</b> |   |
| 00-6200-048-22       | F/M Acetabular Shell, 48mm OD, Cluster    |
| Through ↓            | Through ↓                                 |
| 00-6200-070-22       | F/M Acetabular Shell, 70mm OD, Cluster    |
| 65-6200-048-22       | HA/TCP Acetabular Shell, 48mm OD, Cluster |
| Through ↓            | Through ↓                                 |
| 65-6200-070-22       | HA/TCP Acetabular Shell, 70mm OD, Cluster |
| <b>Multi-Holed</b>   |   |
| 00-6200-036-20       | F/M Acetabular Shell, 36mm OD, Multi      |
| Through ↓            | Through ↓                                 |
| 00-6200-080-20       | F/M Acetabular Shell, 80mm OD, Multi      |
| 65-6200-036-20       | HA/TCP Acetabular Shell, 36mm OD, Multi   |
| Through ↓            | Through ↓                                 |
| 65-6200-080-20       | HA/TCP Acetabular Shell, 80mm OD, Multi   |
| <b>Uni-Holed</b>     |   |
| 00-6200-040-24       | F/M Acetabular Shell, 40mm OD, Uni        |
| Through ↓            | Through ↓                                 |
| 00-6200-070-24       | F/M Acetabular Shell, 70mm OD, Uni        |

### Bone Screws

SCREWS are available in 5mm LENGTH increments  
up to 40mm 6250-SIZE-LENGTH

| Cat. No.       | Description                        |
|----------------|------------------------------------|
| 00-6250-045-15 | Bone Screw, 4.5x15mm, Self-Tapping |
| Through ↓      | Through ↓                          |
| 00-6250-045-40 | Bone Screw, 4.5x40mm, Self-Tapping |
| 00-6250-045-50 | Bone Screw, 4.5x50mm, Self-Tapping |
| 00-6250-045-60 | Bone Screw, 4.5x60mm, Self-Tapping |
| 00-6250-065-15 | Bone Screw, 6.5x15mm, Self-Tapping |
| Through ↓      | Through ↓                          |
| 00-6250-065-40 | Bone Screw, 6.5x40mm, Self-Tapping |
| 00-6250-065-50 | Bone Screw, 6.5x50mm, Self-Tapping |
| 00-6250-065-60 | Bone Screw, 6.5x60mm, Self-Tapping |
| 00-6250-065-70 | Bone Screw, 6.5x70mm, Self-Tapping |
| 00-6250-065-80 | Bone Screw, 6.5x80mm, Self-Tapping |

**Warning:** This device is not approved for screw attachment or fixation to the posterior elements (pedicles) of the cervical, thoracic, or lumbar spine.

### Polyethylene Liner Thickness

| Shell OD (mm) | Polyethylene Liner |      |      |      |      |      |
|---------------|--------------------|------|------|------|------|------|
|               | 22mm               | 26mm | 28mm | 32mm | 36mm | 40mm |
| 36            | 4.3                | N/A  | N/A  | N/A  | N/A  | N/A  |
| 38            | 5.3                | N/A  | N/A  | N/A  | N/A  | N/A  |
| 40            | 6.1                | N/A  | N/A  | N/A  | N/A  | N/A  |
| 42            | 6.1                | 5.3  | N/A  | N/A  | N/A  | N/A  |
| 44            | 7.1                | 6.3  | 5.3  | N/A  | N/A  | N/A  |
| 46            | 8.2                | 6.3  | 6.3  | N/A  | N/A  | N/A  |
| 48            | 9.1                | 7.3  | 6.3  | 5.3  | N/A  | N/A  |
| 50            | 10.1               | 8.3  | 7.3  | 6.3  | 6.7  | N/A  |
| 52            | 10.1               | 8.3  | 7.3  | 6.3  | 6.7  | N/A  |
| 54            | 10.1               | 8.3  | 7.3  | 6.3  | 6.7  | N/A  |
| 56            | 11.3               | 9.4  | 8.4  | 6.4  | 7.8  | N/A  |
| 58            | 12.2               | 10.4 | 9.4  | 7.4  | 8.8  | 6.9  |
| 60            | 13.2               | 11.4 | 10.4 | 8.4  | 9.8  | 7.9  |
| 62            | 14.2               | 12.4 | 11.4 | 9.4  | 10.8 | 8.9  |
| 64            | 15.2               | 13.4 | 12.4 | 10.4 | 11.8 | 9.9  |
| 66            | 16.2               | 14.4 | 13.4 | 11.4 | 12.8 | 10.9 |
| 68            | 17.2               | 15.4 | 14.4 | 12.4 | 13.8 | 11.9 |
| 70            | 18.2               | 16.4 | 15.4 | 13.4 | 14.8 | 12.9 |
| 72            | 19.2               | 17.3 | 16.3 | 14.3 | 15.8 | 13.8 |
| 74            | 20.3               | 18.4 | 17.4 | 15.4 | 16.8 | 14.9 |
| 76            | 21.2               | 19.3 | 18.3 | 16.3 | 17.8 | 15.8 |
| 78            | 22.2               | 20.4 | 19.4 | 17.4 | 18.8 | 16.9 |
| 80            | 23.2               | 21.3 | 20.3 | 18.3 | 19.8 | 17.8 |

### Replacement Liner Locking Rings

| Cat. No. | Description                                |
|----------|--|
| 6201-036 | Modular Cup Replacement Locking Ring, 36mm |
| 6201-038 | Modular Cup Replacement Locking Ring, 38mm |
| 6201-040 | Modular Cup Replacement Locking Ring, 40mm |
| 6201-042 | Modular Cup Replacement Locking Ring, 42mm |
| 6201-044 | Modular Cup Replacement Locking Ring, 44mm |
| 6201-046 | Modular Cup Replacement Locking Ring, 46mm |
| 6201-048 | Modular Cup Replacement Locking Ring, 48mm |
| 6201-050 | Modular Cup Replacement Locking Ring, 50mm |
| 6201-052 | Modular Cup Replacement Locking Ring, 52mm |
| 6201-054 | Modular Cup Replacement Locking Ring, 54mm |
| 6201-056 | Modular Cup Replacement Locking Ring, 56mm |
| 6201-058 | Modular Cup Replacement Locking Ring, 58mm |
| 6201-060 | Modular Cup Replacement Locking Ring, 60mm |
| 6201-062 | Modular Cup Replacement Locking Ring, 62mm |
| 6201-064 | Modular Cup Replacement Locking Ring, 64mm |
| 6201-066 | Modular Cup Replacement Locking Ring, 66mm |
| 6201-068 | Modular Cup Replacement Locking Ring, 68mm |
| 6201-070 | Modular Cup Replacement Locking Ring, 70mm |
| 6201-072 | Modular Cup Replacement Locking Ring, 72mm |
| 6201-074 | Modular Cup Replacement Locking Ring, 74mm |
| 6201-076 | Modular Cup Replacement Locking Ring, 76mm |
| 6201-078 | Modular Cup Replacement Locking Ring, 78mm |
| 6201-080 | Modular Cup Replacement Locking Ring, 80mm |

### Longevity® Highly Crosslinked Polyethylene Liners

Liners are available in 2mm increments  
Shell sizes 50, 52, and 54 use the same size liner.

| Cat. No.       | Description                   |
|----------------|-------------------------------|
| <b>Neutral</b> |                               |
| 00-6305-036-22 | Poly Liner, 36mm OD x 22mm ID |
| Through ↓      | Through ↓                     |
| 00-6305-080-22 | Poly Liner, 80mm OD x 22mm ID |
| 00-6305-042-26 | Poly Liner, 42mm OD x 26mm ID |
| Through ↓      | Through ↓                     |
| 00-6305-080-26 | Poly Liner, 80mm OD x 26mm ID |
| 00-6305-044-28 | Poly Liner, 44mm OD x 28mm ID |
| Through ↓      | Through ↓                     |
| 00-6305-080-28 | Poly Liner, 80mm OD x 28mm ID |
| 00-6305-048-32 | Poly Liner, 48mm OD x 32mm ID |
| Through ↓      | Through ↓                     |
| 00-6305-080-32 | Poly Liner, 80mm OD x 32mm ID |
| 00-6305-050-36 | Poly Liner, 50mm OD x 36mm ID |
| Through ↓      | Through ↓                     |
| 00-6305-080-36 | Poly Liner, 80mm OD x 36mm ID |
| 00-6305-058-40 | Poly Liner, 58mm OD x 40mm ID |
| Through ↓      | Through ↓                     |
| 00-6305-080-40 | Poly Liner, 80mm OD x 40mm ID |

### Longevity Crosslinked Polyethylene Liners

Liners are available in 2mm increments  
Shell sizes 50, 52, and 54 use the same size liner.

| Cat. No.            | Description                               |
|---------------------|---|
| <b>10° Elevated</b> |   |
| 00-6310-036-22      | 10° Elevated Rim Liner, 36mm OD x 22mm ID |
| Through ↓           | Through ↓                                 |
| 00-6310-080-22      | 10° Elevated Rim Liner, 80mm OD x 22mm ID |
| 00-6310-042-26      | 10° Elevated Rim Liner, 42mm OD x 26mm ID |
| Through ↓           | Through ↓                                 |
| 00-6310-080-26      | 10° Elevated Rim Liner, 80mm OD x 26mm ID |
| 00-6310-044-28      | 10° Elevated Rim Liner, 44mm OD x 28mm ID |
| Through ↓           | Through ↓                                 |
| 00-6310-080-28      | 10° Elevated Rim Liner, 80mm OD x 28mm ID |
| 00-6310-048-32      | 10° Elevated Rim Liner, 48mm OD x 32mm ID |
| Through ↓           | Through ↓                                 |
| 00-6310-080-32      | 10° Elevated Rim Liner, 80mm OD x 32mm ID |
| 00-6310-058-36      | 10° Elevated Rim Liner, 58mm OD x 36mm ID |
| Through ↓           | Through ↓                                 |
| 00-6310-080-36      | 10° Elevated Rim Liner, 80mm OD x 36mm ID |
| <b>20° Elevated</b> |   |
| 00-6320-036-22      | 20° Elevated Rim Liner, 36mm OD x 22mm ID |
| Through ↓           | Through ↓                                 |
| 00-6320-080-22      | 20° Elevated Rim Liner, 80mm OD x 22mm ID |
| 00-6320-042-26      | 20° Elevated Rim Liner, 42mm OD x 26mm ID |
| Through ↓           | Through ↓                                 |
| 00-6320-080-26      | 20° Elevated Rim Liner, 80mm OD x 26mm ID |

For more information on the *Trilogy Longevity* Constrained Liners and Instrumentation, refer to Surgical Technique, Lit # 97-6334-002-00.

|                |   |
|----------------|---|
| 00-6320-044-28 | 20° Elevated Rim Liner, 44mm OD x 28mm ID |
| Through ↓      | Through ↓                                 |
| 00-6320-080-28 | 20° Elevated Rim Liner, 80mm OD x 28mm ID |
| 00-6320-048-32 | 20° Elevated Rim Liner, 48mm OD x 32mm ID |
| Through ↓      | Through ↓                                 |
| 00-6320-080-32 | 20° Elevated Rim Liner, 80mm OD x 32mm ID |

#### 7mm Offset

|                |                               |
|----------------|-------------------------------|
| 00-6341-040-22 | 7mm Offset, 40mm OD x 22mm ID |
| Through ↓      | Through ↓                     |
| 00-6341-070-22 | 7mm Offset, 70mm OD x 22mm ID |
| 00-6341-042-26 | 7mm Offset, 42mm OD x 26mm ID |
| Through ↓      | Through ↓                     |
| 00-6341-070-26 | 7mm Offset, 70mm OD x 26mm ID |
| 00-6341-044-28 | 7mm Offset, 44mm OD x 28mm ID |
| Through ↓      | Through ↓                     |
| 00-6341-070-28 | 7mm Offset, 70mm OD x 28mm ID |
| 00-6341-048-32 | 7mm Offset, 48mm OD x 32mm ID |
| Through ↓      | Through ↓                     |
| 00-6341-070-32 | 7mm Offset, 70mm OD x 32mm ID |

## Standard Polyethylene Liners

Liners are available in 2mm increments  
Shell sizes 50, 52, and 54 use the same size liner.

| Cat. No. | Description |
|----------|-------------|
|----------|-------------|

#### Neutral

|                |                               |
|----------------|-------------------------------|
| 00-6105-036-22 | Poly Liner, 36mm OD x 22mm ID |
| Through ↓      | Through ↓                     |
| 00-6105-080-22 | Poly Liner, 80mm OD x 22mm ID |
| 00-6105-042-26 | Poly Liner, 42mm OD x 26mm ID |
| Through ↓      | Through ↓                     |
| 00-6105-080-26 | Poly Liner, 80mm OD x 26mm ID |
| 00-6105-044-28 | Poly Liner, 44mm OD x 28mm ID |
| Through ↓      | Through ↓                     |
| 00-6105-080-28 | Poly Liner, 80mm OD x 28mm ID |
| 00-6105-048-32 | Poly Liner, 48mm OD x 32mm ID |
| Through ↓      | Through ↓                     |
| 00-6105-080-32 | Poly Liner, 80mm OD x 32mm ID |

#### 10° Elevated

|                |   |
|----------------|---|
| 00-6110-036-22 | 10° Elevated Rim Liner, 36mm OD x 22mm ID |
| Through ↓      | Through ↓                                 |
| 00-6110-080-22 | 10° Elevated Rim Liner, 80mm OD x 22mm ID |
| 00-6110-042-26 | 10° Elevated Rim Liner, 42mm OD x 26mm ID |
| Through ↓      | Through ↓                                 |
| 00-6110-080-26 | 10° Elevated Rim Liner, 80mm OD x 26mm ID |
| 00-6110-044-28 | 10° Elevated Rim Liner, 44mm OD x 28mm ID |
| Through ↓      | Through ↓                                 |
| 00-6110-080-28 | 10° Elevated Rim Liner, 80mm OD x 28mm ID |
| 00-6110-048-32 | 10° Elevated Rim Liner, 48mm OD x 32mm ID |
| Through ↓      | Through ↓                                 |
| 00-6110-080-32 | 10° Elevated Rim Liner, 80mm OD x 32mm ID |

#### 20° Elevated

|                |   |
|----------------|---|
| 00-6120-036-22 | 20° Elevated Rim Liner, 36mm OD x 22mm ID |
| Through ↓      | Through ↓                                 |
| 00-6120-080-22 | 20° Elevated Rim Liner, 80mm OD x 22mm ID |
| 00-6120-042-26 | 20° Elevated Rim Liner, 42mm OD x 26mm ID |
| Through ↓      | Through ↓                                 |
| 00-6120-080-26 | 20° Elevated Rim Liner, 80mm OD x 26mm ID |
| 00-6120-044-28 | 20° Elevated Rim Liner, 44mm OD x 28mm ID |
| Through ↓      | Through ↓                                 |
| 00-6120-080-28 | 20° Elevated Rim Liner, 80mm OD x 28mm ID |
| 00-6120-048-32 | 20° Elevated Rim Liner, 48mm OD x 32mm ID |
| Through ↓      | Through ↓                                 |
| 00-6120-080-32 | 20° Elevated Rim Liner, 80mm OD x 32mm ID |

## Standard Polyethylene Liners

Liners are available in 2mm increments  
Shell sizes 50, 52, and 54 use the same size liner.

| Cat. No. | Description |
|----------|-------------|
|----------|-------------|

#### 7mm Offset

|                |                               |
|----------------|-------------------------------|
| 00-6141-040-22 | 7mm Offset, 40mm OD x 22mm ID |
| Through ↓      | Through ↓                     |
| 00-6141-070-22 | 7mm Offset, 70mm OD x 22mm ID |
| 00-6141-042-26 | 7mm Offset, 42mm OD x 26mm ID |
| Through ↓      | Through ↓                     |
| 00-6141-070-26 | 7mm Offset, 70mm OD x 26mm ID |
| 00-6141-044-28 | 7mm Offset, 44mm OD x 28mm ID |
| Through ↓      | Through ↓                     |
| 00-6141-070-28 | 7mm Offset, 70mm OD x 28mm ID |
| 00-6141-048-32 | 7mm Offset, 48mm OD x 32mm ID |
| Through ↓      | Through ↓                     |
| 00-6141-070-32 | 7mm Offset, 70mm OD x 32mm ID |

#### Eccentric

|                |   |
|----------------|---|
| 00-6151-000-28 | Eccentric Liner Set<br>(Includes all 28mm ID liners listed below) |
| 00-6151-064-28 | Eccentric Liner, 64mm OD x 28mm ID                                |
| Through ↓      | Through ↓   |
| 00-6151-080-28 | Eccentric Liner, 80mm OD x 28mm ID                                |

|                |   |
|----------------|---|
| 00-6151-000-32 | Eccentric Liner Set<br>(Includes all 32mm ID liners listed below) |
| 00-6151-064-32 | Eccentric Liner, 64mm OD x 32mm ID                                |
| Through ↓      | Through ↓   |
| 00-6151-080-32 | Eccentric Liner, 80mm OD x 32mm ID                                |

#### Oblique

|                |  |
|----------------|--|
| 00-6152-000-01 | Oblique Liner Set<br>(Includes all 22mm liners listed below) |
| 00-6152-038-22 | Oblique Liner, 38mm OD x 22mm ID                             |
| Through ↓      | Through ↓  |
| 00-6152-048-22 | Oblique Liner, 48mm OD x 22mm ID                             |

|                |  |
|----------------|--|
| 00-6152-000-02 | Oblique Liner Set<br>(Includes all 28mm liners listed below) |
|----------------|--|

|                |                                  |
|----------------|----------------------------------|
| 00-6152-050-28 | Oblique Liner, 50mm OD x 28mm ID |
| Through ↓      | Through ↓                        |
| 00-6152-062-28 | Oblique Liner, 62mm OD x 28mm ID |

|                |  |
|----------------|--|
| 00-6152-000-03 | Oblique Liner Set<br>(Includes all 32mm liners listed below) |
| 00-6152-056-32 | Oblique Liner, 56mm OD x 32mm ID                             |
| Through ↓      | Through ↓  |
| 00-6152-062-32 | Oblique Liner, 62mm OD x 32mm ID                             |

## Trilogy Longevity Constrained Liners (Includes Constraining Ring)

| Cat. No.       | Description   |
|----------------|---|
| 00-6334-050-28 | <i>Trilogy Longevity</i> Constrained Liner<br>50/52/54mm OD x 28mm ID |
| 00-6334-056-32 | <i>Trilogy Longevity</i> Constrained Liner<br>56mm OD x 32mm ID       |
| 00-6334-058-32 | <i>Trilogy Longevity</i> Constrained Liner<br>58mm OD x 32mm ID       |
| 00-6334-060-36 | <i>Trilogy Longevity</i> Constrained Liner<br>60mm OD x 36mm ID       |
| 00-6334-062-36 | <i>Trilogy Longevity</i> Constrained Liner<br>62mm OD x 36mm ID       |
| 00-6334-064-36 | <i>Trilogy Longevity</i> Constrained Liner<br>64mm OD x 36mm ID       |
| 00-6334-066-36 | <i>Trilogy Longevity</i> Constrained Liner<br>66mm OD x 36mm ID       |
| 00-6334-068-36 | <i>Trilogy Longevity</i> Constrained Liner<br>68mm OD x 36mm ID       |
| 00-6334-070-36 | <i>Trilogy Longevity</i> Constrained Liner<br>70mm OD x 36mm ID       |
| 00-6334-072-36 | <i>Trilogy Longevity</i> Constrained Liner<br>72mm OD x 36mm ID       |
| 00-6334-074-36 | <i>Trilogy Longevity</i> Constrained Liner<br>74mm OD x 36mm ID       |
| 00-6334-076-36 | <i>Trilogy Longevity</i> Constrained Liner<br>76mm OD x 36mm ID       |
| 00-6334-078-36 | <i>Trilogy Longevity</i> Constrained Liner<br>78mm OD x 36mm ID       |
| 00-6334-080-36 | <i>Trilogy Longevity</i> Constrained Liner<br>80mm OD x 36mm ID       |

### Trilogy Longevity Constrained Liner Instruments

| Cat. No.       | Description   |
|----------------|---|
| 00-6144-000-00 | <i>Trilogy Longevity</i> Constrained Liner Instrument Kit (includes one each of the items listed below) |
| 00-6144-095-00 | <i>Trilogy Longevity</i> Constrained Liner Instrument Case (includes base, tray, and lid)               |
| 00-6144-050-28 | <i>Trilogy Longevity</i> Constrained Liner Provisional OD x ID 50, 52, 54 x 28                          |
| 00-6144-056-32 | <i>Trilogy Longevity</i> Constrained Liner Provisional OD x ID 56 x 32                                  |
| 00-6144-058-32 | <i>Trilogy Longevity</i> Constrained Liner Provisional OD x ID 58 x 32                                  |
| 00-6144-060-36 | <i>Trilogy Longevity</i> Constrained Liner Provisional OD x ID 60 x 36                                  |
| 00-6144-062-36 | <i>Trilogy Longevity</i> Constrained Liner Provisional OD x ID 62 x 36                                  |
| 00-6144-064-36 | <i>Trilogy Longevity</i> Constrained Liner Provisional OD x ID 64 x 36                                  |
| 00-6144-066-36 | <i>Trilogy Longevity</i> Constrained Liner Provisional OD x ID 66 x 36                                  |
| 00-6144-068-36 | <i>Trilogy Longevity</i> Constrained Liner Provisional OD x ID 68 x 36                                  |
| 00-6144-070-36 | <i>Trilogy Longevity</i> Constrained Liner Provisional OD x ID 70 x 36                                  |
| 9340-00-000    | Impactor Handle   |
| 00-6144-001-28 | Constraining Ring Impactor 28mm   |
| 00-6144-001-32 | Constraining Ring Impactor 32mm   |
| 00-6144-001-36 | Constraining Ring Impactor 36mm   |
| 00-6144-002-28 | Liner Impactor Head 28mm  |
| 00-6144-002-32 | Liner Impactor Head 32mm  |
| 00-6144-002-36 | Liner Impactor Head 36mm  |

### Trilogy Longevity Constrained Liner Jumbo Liner Provisionals

| Cat. No.       | Description  |
|----------------|--|
| 00-6144-000-01 | <i>Trilogy Longevity</i> Constrained Liner Jumbo Provisionals Kit (includes one each of the items listed below) (fits within case #00-6144-095-00) |
| 00-6144-072-36 | <i>Trilogy Longevity</i> Constrained Liner Provisional OD x ID 72 x 36   |
| 00-6144-074-36 | <i>Trilogy Longevity</i> Constrained Liner Provisional OD x ID 74 x 36   |
| 00-6144-076-36 | <i>Trilogy Longevity</i> Constrained Liner Provisional OD x ID 76 x 36   |
| 00-6144-078-36 | <i>Trilogy Longevity</i> Constrained Liner Provisional OD x ID 78 x 36   |
| 00-6144-080-36 | <i>Trilogy Longevity</i> Constrained Liner Provisional OD x ID 80 x 36   |

### Trilogy Bone Screws

| Cat. No.       | Description                       |
|----------------|-----------------------------------|
| 00-6250-045-15 | Bone Screw 4.5mm x 15mm, Self-Tap |
| 00-6250-045-20 | Bone Screw 4.5mm x 20mm, Self-Tap |
| 00-6250-045-25 | Bone Screw 4.5mm x 25mm, Self-Tap |
| 00-6250-045-30 | Bone Screw 4.5mm x 30mm, Self-Tap |
| 00-6250-045-35 | Bone Screw 4.5mm x 35mm, Self-Tap |
| 00-6250-045-40 | Bone Screw 4.5mm x 40mm, Self-Tap |
| 00-6250-045-50 | Bone Screw 4.5mm x 50mm, Self-Tap |
| 00-6250-045-60 | Bone Screw 4.5mm x 60mm, Self-Tap |
| 00-6250-065-15 | Bone Screw 6.5mm x 15mm, Self-Tap |
| 00-6250-065-20 | Bone Screw 6.5mm x 20mm, Self-Tap |
| 00-6250-065-25 | Bone Screw 6.5mm x 25mm, Self-Tap |
| 00-6250-065-30 | Bone Screw 6.5mm x 30mm, Self-Tap |
| 00-6250-065-35 | Bone Screw 6.5mm x 35mm, Self-Tap |
| 00-6250-065-40 | Bone Screw 6.5mm x 40mm, Self-Tap |
| 00-6250-065-50 | Bone Screw 6.5mm x 50mm, Self-Tap |
| 00-6250-065-60 | Bone Screw 6.5mm x 60mm, Self-Tap |

Warning: This device is not approved for screw attachment or fixation to the posterior elements (pedicles) of the cervical, thoracic, or lumbar spine.

### Trilogy Holed Instruments

| Cat. No.       | Description   |
|----------------|---|
| 00-6260-099-02 | <i>Trilogy</i> Holed Instrument Set (Set includes one each of the items listed below) |
| 00-6260-085-01 | Instrument Case (includes base and lid)   |
| 00-6260-002-00 | Flex Shaft with Modular Connector   |
| 00-6260-003-01 | Drill Bit, 15mm   |
| 00-6260-003-02 | Drill Bit, 30mm   |
| 00-6260-003-03 | Drill Bit, 45mm   |
| 00-6260-006-00 | Drill Guide   |
| 00-6260-007-01 | Tap, 4.5mm  |
| 00-6260-007-02 | Tap, 6.5mm  |
| 00-6260-008-01 | Tap Guide, 4.5mm  |
| 00-6260-008-02 | Tap Guide, 6.5mm  |
| 00-6260-10-00  | Tap Handle  |
| 00-6260-013-00 | Screw Holding Forceps, 15°  |
| 00-6260-014-00 | Screw Holding Forceps, 45°  |
| 00-6260-024-00 | Straight Screwdriver  |
| 00-6260-025-00 | Universal Screwdriver   |
| 00-6260-026-00 | Modular Universal Handle  |
| 00-6611-098-00 | Depth Gauge   |
| 00-4215-200-00 | Screwdriver Torque Limiter (ordered separately, not in kit)                           |

### Liner Removal Instruments

(Included in *Trilogy Non-Holed Instrument Set*, 00-6260-099-01)

| Cat. No.       | Description        |
|----------------|--------------------|
| 00-6260-030-01 | Disassembly Device |
| 00-6260-035-01 | Liner Extractor    |
| 00-6260-040-00 | Liner Elevator     |

### Liner Provisionals

Liner Provisionals are available in 2mm increments. Shell sizes 50, 52, and 54 use the same size provisional liner. Locking Screws are ordered separately.

#### Polar Screws (Optional)

| Cat. No.       | Description   |
|----------------|---|
| 00-6270-099-10 | Polar Screw Set (Includes seven each of item below) |
| 00-6270-001-00 | Polar Screw   |

#### Liner Provisionals 22mm

| Cat. No.       | Description   |
|----------------|---|
| 00-6270-099-22 | Neutral Provisional Liner Set - 22mm (Includes one each of items below)         |
| 00-6270-040-22 | Neutral Provisional Liner, 22mm   |
| Through ↓      | Through ↓   |
| 00-6270-070-22 | Neutral Provisional Liner, 22mm   |
| 00-6271-099-22 | 10° Elev. Provisional Liner Set - 22mm (Includes one each of items below)       |
| 00-6271-040-22 | 10° Elev. Provisional Liner, 22mm   |
| Through ↓      | Through ↓   |
| 00-6271-070-22 | 10° Elev. Provisional Liner, 22mm   |
| 00-6272-099-22 | 20° Elev. Provisional Liner Set - 22mm (Includes one each of items below)       |
| 00-6272-040-22 | 20° Elev. Provisional Liner, 22mm   |
| Through ↓      | Through ↓   |
| 00-6272-070-22 | 20° Elev. Provisional Liner, 22mm   |
| 00-6270-199-22 | Neutral Jumbo Provisional Liner Set - 22mm (Includes one each of items below)   |
| 00-6270-072-22 | Neutral Jumbo Provisional Liner, 22mm   |
| Through ↓      | Through ↓   |
| 00-6270-080-22 | Neutral Jumbo Provisional Liner, 22mm   |
| 00-6271-199-22 | 10° Elev. Jumbo Provisional Liner Set - 22mm (Includes one each of items below) |
| 00-6271-072-22 | 10° Elev. Jumbo Provisional Liner, 22mm   |
| Through ↓      | Through ↓   |
| 00-6271-080-22 | 10° Jumbo Provisional Liner, 22mm   |
| 00-6272-199-22 | 20° Elev. Jumbo Provisional Liner Set - 22mm (Includes one each of items below) |
| 00-6272-072-22 | 20° Elev. Provisional Liner, 22mm   |
| Through ↓      | Through ↓   |
| 00-6272-080-22 | 20° Elev. Provisional Liner, 22mm   |

**Liner Provisionals 26mm**

| Cat. No.       | Description  |
|----------------|--|
| 00-6270-099-26 | Neutral Provisional Liner Set - 26mm<br>(Includes one each of items below)         |
| 00-6270-042-26 | Neutral Provisional Liner, 26mm  |
| Through ↓      | Through ↓  |
| 00-6270-070-26 | Neutral Provisional Liner, 26mm  |
| 00-6271-099-26 | 10° Elev. Provisional Liner Set - 26mm<br>(Includes one each of items below)       |
| 00-6271-042-26 | 10° Elev. Provisional Liner, 26mm  |
| Through ↓      | Through ↓  |
| 00-6271-070-26 | 10° Elev. Provisional Liner, 26mm  |
| 00-6272-099-26 | 20° Elev. Provisional Liner Set - 26mm<br>(Includes one each of items below)       |
| 00-6272-042-26 | 20° Elev. Provisional Liner, 26mm  |
| Through ↓      | Through ↓  |
| 00-6272-070-26 | 20° Elev. Provisional Liner, 26mm  |
| 00-6270-199-26 | Neutral Jumbo Provisional Liner Set - 26mm<br>(Includes one each of items below)   |
| 00-6270-072-26 | Neutral Jumbo Provisional Liner, 26mm  |
| Through ↓      | Through ↓  |
| 00-6270-080-26 | Neutral Jumbo Provisional Liner, 26mm  |
| 00-6271-199-26 | 10° Elev. Jumbo Provisional Liner Set - 26mm<br>(Includes one each of items below) |
| 00-6271-072-26 | 10° Elev. Jumbo Provisional Liner, 26mm  |
| Through ↓      | Through ↓  |
| 00-6271-080-26 | 10° Elev. Jumbo Provisional Liner, 26mm  |
| 00-6272-199-26 | 20° Elev. Jumbo Provisional Liner Set - 26mm<br>(Includes one each of items below) |
| 00-6272-072-26 | 20° Elev. Jumbo Provisional Liner, 26mm  |
| Through ↓      | Through ↓  |
| 00-6272-080-26 | 20° Elev. Jumbo Provisional Liner, 26mm  |

**Liner Provisionals 28mm**

| Cat. No.       | Description  |
|----------------|--|
| 00-6270-099-28 | Neutral Provisional Liner Set - 28mm<br>(Includes one each of items below)       |
| 00-6270-044-28 | Neutral Provisional Liner, 28mm  |
| Through ↓      | Through ↓  |
| 00-6270-070-28 | Neutral Provisional Liner, 28mm  |
| 00-6271-099-28 | 10° Elev. Provisional Liner Set - 28mm<br>(Includes one each of items below)     |
| 00-6271-044-28 | 10° Elev. Provisional Liner, 28mm  |
| Through ↓      | Through ↓  |
| 00-6271-070-28 | 10° Elev. Provisional Liner, 28mm  |
| 00-6272-099-28 | 20° Elev. Provisional Liner Set - 28mm<br>(Includes one each of items below)     |
| 00-6272-044-28 | 20° Elev. Provisional Liner, 28mm  |
| Through ↓      | Through ↓  |
| 00-6272-070-28 | 20° Elev. Provisional Liner, 28mm  |
| 00-6270-199-28 | Neutral Jumbo Provisional Liner Set - 28mm<br>(Includes one each of items below) |
| 00-6270-072-28 | Neutral Jumbo Provisional Liner, 28mm  |
| Through ↓      | Through ↓  |
| 00-6270-080-28 | Neutral Jumbo Provisional Liner, 28mm  |

**Liner Provisionals 28mm (cont.)**

| Cat. No.       | Description  |
|----------------|--|
| 00-6271-199-28 | 10° Elev. Jumbo Provisional Liner Set - 28mm<br>(Includes one each of items below) |
| 00-6271-072-28 | 10° Elev. Jumbo Provisional Liner, 28mm  |
| Through ↓      | Through ↓  |
| 00-6271-080-28 | 10° Elev. Jumbo Provisional Liner, 28mm  |
| 00-6272-199-28 | 20° Elev. Jumbo Provisional Liner Set - 28mm<br>(Includes one each of items below) |
| 00-6272-072-28 | 20° Elev. Jumbo Provisional Liner, 28mm  |
| Through ↓      | Through ↓  |
| 00-6272-080-28 | 20° Elev. Jumbo Provisional Liner, 28mm  |

**Liner Provisionals 32mm**

| Cat. No.       | Description  |
|----------------|--|
| 00-6270-099-32 | Neutral Provisional Liner Set - 32mm<br>(Includes one each of items below)         |
| 00-6270-048-32 | Neutral Provisional Liner, 32mm  |
| Through ↓      | Through ↓  |
| 00-6270-070-32 | Neutral Provisional Liner, 32mm  |
| 00-6271-099-32 | 10° Elev. Provisional Liner Set - 32mm<br>(Includes one each of items below)       |
| 00-6271-048-32 | 10° Elev. Provisional Liner, 32mm  |
| Through ↓      | Through ↓  |
| 00-6271-070-32 | 10° Elev. Provisional Liner, 32mm  |
| 00-6272-099-32 | 20° Provisional Liner Set - 32mm<br>(Includes one each of items below)             |
| 00-6272-048-32 | 20° Elev. Provisional Liner, 32mm  |
| Through ↓      | Through ↓  |
| 00-6272-070-32 | 20° Elev. Provisional Liner, 32mm  |
| 00-6270-199-32 | Neutral Jumbo Provisional Liner Set - 32mm<br>(Includes one each of items below)   |
| 00-6270-072-32 | Neutral Jumbo Provisional Liner, 32mm  |
| Through ↓      | Through ↓  |
| 00-6270-080-32 | Neutral Jumbo Provisional Liner, 32mm  |
| 00-6271-199-32 | 10° Elev. Jumbo Provisional Liner Set - 32mm<br>(Includes one each of items below) |
| 00-6271-072-32 | 10° Elev. Jumbo Provisional Liner, 32mm  |
| Through ↓      | Through ↓  |
| 00-6271-080-32 | 10° Elev. Jumbo Provisional Liner, 32mm  |
| 00-6272-199-32 | 20° Elev. Jumbo Provisional Liner Set - 32mm<br>(Includes one each of items below) |
| 00-6272-072-32 | 20° Elev. Jumbo Provisional Liner, 32mm  |
| Through ↓      | Through ↓  |
| 00-6272-080-32 | 20° Elev. Jumbo Provisional Liner, 32mm  |

**Liner Provisionals/Instruments 36mm**

| Cat. No.       | Description  |
|----------------|--|
| 00-6270-099-36 | Neutral Provisional Liner Set - 36mm<br>(Includes one each of items below)   |
| 00-6261-090-00 | Case (Includes Base and Lid)   |
| 00-6270-050-36 | Neutral Provisional Liner, 36mm  |
| Through ↓      | Through ↓  |
| 00-6270-080-36 | Neutral Provisional Liner, 36mm  |
| 00-7695-036-01 | 6° Femoral Head Provisional  |
| Through ↓      | Through ↓  |
| 00-7695-036-05 | 6° Femoral Head Provisional  |
| 00-7895-036-01 | 12/14 Femoral Head Provisional   |
| Through ↓      | Through ↓  |
| 00-7895-036-05 | 12/14 Femoral Head Provisional   |
| 00-6260-017-36 | Liner Insertter  |
| 00-6271-099-36 | 10° Elev. Provisional Liner Set - 36mm<br>(Includes one each of items below) |
| 00-6261-090-00 | Case (Includes Base and Lid)   |
| 00-6271-058-36 | 10° Elev. Provisional Liner, 36mm  |
| Through ↓      | Through ↓  |
| 00-6271-080-36 | 10° Elev. Provisional Liner, 36mm  |
| 00-7695-036-01 | 6° Femoral Head Provisional  |
| Through ↓      | Through ↓  |
| 00-7695-036-05 | 6° Femoral Head Provisional  |
| 00-7895-036-01 | 12/14 Femoral Head Provisional   |
| Through ↓      | Through ↓  |
| 00-7895-036-05 | 12/14 Femoral Head Provisional   |
| 00-6260-017-36 | Liner Insertter  |

**Liner Provisionals/Instruments 40mm**

| Cat. No.       | Description   |
|----------------|---|
| 00-6270-099-40 | Neutral Instrument Set - 40mm<br>(Includes one each of items below) |
| 00-6270-002-00 | Case (Includes Base and Lid)  |
| 00-6270-058-40 | Neutral Provisional Liner, 40mm                                     |
| Through ↓      | Through ↓   |
| 00-6270-080-40 | Neutral Provisional Liner, 40mm                                     |
| 00-7895-040-01 | 12/14 Femoral Head Provisional                                      |
| Through ↓      | Through ↓   |
| 00-7895-040-05 | 12/14 Femoral Head Provisional                                      |
| 00-6260-017-40 | Liner Insertter   |



### Liner Provisionals 7mm Offset

| Cat. No.                    | Description   |
|-----------------------------|---|
| 00-6275-099-22              | 7mm Offset Provisional Liner Set - 22mm<br>(Includes one each of items below) |
| 00-6275-040-22<br>Through ↓ | 7mm Offset Provisional Liner, 22mm<br>Through ↓                               |
| 00-6275-070-22              | 7mm Offset Provisional Liner, 22mm  |
| 00-6275-099-26              | 7mm Offset Provisional Liner Set - 26mm<br>(Includes one each of items below) |
| 00-6275-042-26<br>Through ↓ | 7mm Offset Provisional Liner, 26mm<br>Through ↓                               |
| 00-6275-070-26              | 7mm Offset Provisional Liner, 26mm  |
| 00-6275-099-28              | 7mm Offset Provisional Liner Set - 28mm<br>(Includes one each of items below) |
| 00-6275-044-28<br>Through ↓ | 7mm Offset Provisional Liner, 28mm<br>Through ↓                               |
| 00-6275-070-28              | 7mm Offset Provisional Liner, 28mm  |
| 00-6275-099-32              | 7mm Offset Provisional Liner Set - 32mm<br>(Includes one each of items below) |
| 00-6275-048-32<br>Through ↓ | 7mm Offset Provisional Liner, 32mm<br>Through ↓                               |
| 00-6275-070-32              | 7mm Offset Provisional Liner, 32mm  |

### Liner Provisionals - Eccentric

| Cat. No.                    | Description  |
|-----------------------------|--|
| 00-6171-099-28              | Eccentric Provisional Liner Set - 28mm<br>(Includes one each of items below) |
| 00-6171-064-28<br>Through ↓ | Eccentric Provisional Liner, 28mm<br>Through ↓                               |
| 00-6161-080-28              | Eccentric Provisional Liner, 28mm  |
| 00-6171-099-32              | Eccentric Provisional Liner Set - 32mm<br>(Includes one each of items below) |
| 00-6171-064-32<br>Through ↓ | Eccentric Provisional Liner, 32mm<br>Through ↓                               |
| 00-6171-080-32              | Eccentric Provisional Liner, 32mm  |
| 00-6260-019-00              | Eccentric Liner Inserter Handle  |
| 00-6260-019-28              | Eccentric Liner Inserter Head, 28mm  |
| 00-6260-019-32              | Eccentric Liner Inserter Head, 32mm  |

### Liner Provisionals - Oblique

| Cat. No.                    | Description  |
|-----------------------------|--|
| 00-6172-099-22              | Oblique Provisional Liner Set - 22mm<br>(Includes one each of items below) |
| 00-6172-038-22<br>Through ↓ | Oblique Provisional Liner, 22mm<br>Through ↓                               |
| 00-6162-048-22              | Oblique Provisional Liner, 22mm  |
| 00-6172-099-28              | Oblique Provisional Liner Set - 28mm<br>(Includes one each of items below) |
| 00-6172-050-28<br>Through ↓ | Oblique Provisional Liner, 28mm<br>Through ↓                               |
| 00-6172-062-28              | Oblique Provisional Liner, 28mm  |
| 00-6172-099-32              | Oblique Provisional Liner Set - 32mm<br>(Includes one each of items below) |
| 00-6172-056-32<br>Through ↓ | Oblique Provisional Liner, 32mm<br>Through ↓                               |
| 00-6172-062-32              | Oblique Provisional Liner, 32mm  |

### Liner Removal Instruments

(Included in *Trilogy* Non-Holed Instrument Set,  
**00-6260-099-01**)

| Cat. No.       | Description        |
|----------------|--------------------|
| 00-6260-030-01 | Disassembly Device |
| 00-6260-035-01 | Liner Extractor    |
| 00-6260-040-00 | Liner Elevator     |

### Modular Acetabular Instrument Set

| Cat. No.       | Description  |
|----------------|--|
| 00-6270-099-01 | Instrument Set<br>(Includes one each of items below) |
| 00-6275-010-00 | Case (Includes Base and Lid)                         |
| 00-6260-030-01 | Disassembly Device                                   |
| 00-6260-035-01 | Liner Extractor                                      |
| 00-6260-040-00 | Liner Elevator                                       |
| 00-6260-018-00 | Cup Positioner                                       |
| 00-6260-047-00 | Gunsight Alignment Guide                             |
| 00-5785-079-00 | Tibial Alignment Rod                                 |
| 00-6260-002-00 | Flex Shaft w/Modular Connector                       |
| 00-6260-003-01 | Drill Bit, 15mm                                      |
| 00-6260-003-02 | Drill Bit, 30mm                                      |
| 00-6260-003-03 | Drill Bit, 45mm                                      |
| 00-6260-006-00 | Drill Guide  |
| 00-6260-007-01 | Tap, 4.5mm   |
| 00-6260-007-02 | Tap, 6.5mm   |
| 00-6260-008-01 | Tap Guide, 4.5mm                                     |
| 00-6260-008-02 | Tap Guide, 6.5mm                                     |
| 00-6260-010-00 | Tap Handle   |
| 00-6260-013-00 | Screw Holding Forceps, 15"                           |
| 00-6260-014-00 | Screw Holding Forceps, 45"                           |
| 00-6260-024-00 | Straight Screwdriver                                 |
| 00-6260-025-00 | Universal Screwdriver                                |
| 00-6260-026-00 | Modular Universal Handle                             |
| 00-6611-098-00 | Depth Gauge  |
| 00-4215-200-00 | Torque Limiter                                       |

### Optional/Additional Instruments

| Cat. No.       | Description                         |
|----------------|-------------------------------------|
| 00-6260-046-00 | A-Frame Alignment Guide             |
| 00-6611-003-01 | Drill Bit, Short-Flexible (HGP II)  |
| 00-6611-003-02 | Drill Bit, Medium-Flexible (HGP II) |
| 00-6611-003-03 | Drill Bit, Long-Flexible (HGP II)   |
| 00-6260-015-11 | Replaceable Poly Cap (6 ea.)        |
| 00-6260-015-12 | Replaceable Metal Cap (1 ea.)       |
| 00-6260-016-22 | Liner Inserter (Vacuum)             |
| 00-6260-016-26 | Liner Inserter (Vacuum)             |
| 00-6260-016-28 | Liner Inserter (Vacuum)             |
| 00-6260-016-32 | Liner Inserter (Vacuum)             |
| 00-6260-017-22 | Liner Inserter (Mechanical)         |
| 00-6260-017-26 | Liner Inserter (Mechanical)         |
| 00-6260-017-28 | Liner Inserter (Mechanical)         |
| 00-6260-017-32 | Liner Inserter (Mechanical)         |
| 00-6260-017-36 | Liner Inserter (Mechanical)         |
| 00-6260-017-40 | Liner Inserter (Mechanical)         |
| 00-6611-010-00 | HGP II Cup Pusher                   |
| 00-6275-019-00 | Rim Impactor Tray                   |
| 00-6275-025-00 | Jumbo Provisional Half Case         |

### Acetabular Reamers

| Cat. No.       | Description  |
|----------------|--|
| 00-1206-090-01 | Acetabular Reamer System<br>Standard Kit (36 - 70mm) |
| 00-1206-090-02 | Acetabular Reamer System<br>Jumbo Kit (71 - 80mm)    |

### Rim Impactors

Available in 4mm increments.

| Cat. No.                    | Description  |
|-----------------------------|--|
| 00-6270-099-09              | Rim Impactor Set<br>(Includes one each of items below) |
| 00-6270-020-00              | Case (Includes Base and Lid)                           |
| 00-6270-038-00<br>Through ↓ | Rim Impactor, 38/40mm<br>Through ↓                     |
| 00-6270-078-00              | Rim Impactor, 78/80mm                                  |

### Shell Provisionals

| Cat. No.                    | Description   |
|-----------------------------|---|
| 00-6270-099-02              | Shell Provisional Set (Includes one each<br>of items below)       |
| 00-6275-015-00              | Case (Includes Base and Lid)                                      |
| 00-6242-040-00<br>Through ↓ | Shell Provisional, 40mm<br>Through ↓                              |
| 00-6242-070-00              | Shell Provisional, 70mm   |
| 00-6270-199-02              | Jumbo Shell Provisional Set<br>(Includes one each of items below) |
| 00-6275-018-00              | Tray  |
| 00-6242-072-00<br>Through ↓ | Shell Provisional, 72mm<br>Through ↓                              |
| 00-6242-080-00              | Shell Provisional, 80mm   |



Please refer to package insert for complete product information, including contraindications, warnings, precautions, and adverse effects.

Contact your Zimmer Representative or visit us at [www.zimmer.com](http://www.zimmer.com)

